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
## The unparalleled accuracy of hyperspectral analysis

Spectroscopy in near infrared basing on hyperspectral imaging as an instrument to promptly evaluate antioxidant capacity of red wines

Consumers' taste and preferences about wine have changed over the years as a response to new fine taste and possible health issues.

The fact that increased consumption of products rich in natural polyphenols linking to a risk of cardiovascular disease has drawn particular attention to red wines and their relevant content of biologically active compounds. Among them, phenolic compounds impact the most on wine organoleptic characteristics and deliver great benefits in terms of cardiovascular disease. Consumers' desire for quality products that also offer benefits to their health and the growing

number of new wines on the market, have triggered the need to develop increasingly accurate analytical instruments, which are fast and useful to determine the quantity of those bioactive compounds. To this day, there are different analytical techniques we can use to analyse the chemical composition and quality of wine products. In particular, UV-Vis-NIR spectroscopy stands out as one of the best to carry out routine analysis thanks to its fast analysis and low cost.

From a strictly industrial point of view, spectrometric analysis can be carried out by using hyperspectral systems that combine the spectroscopical potential with vision systems thus acquiring a spectrum per every pixel of the image.

Hyperspectral technology operates with both visible light, in spectral band 400-1000nm, and infrared light

in spectral band 900-1700nm. In both cases, light absorption induced by the vibration of the molecules is studied and provides relevant information on the chemistry of the samples analyzed.

Hyperspectral cameras SPECIM FX10, with spectral band 40-1000nm, and FX17 operating with spectral band 900-1700nm, have been used for wine. During the measurements, the concentration of resveratrol in samples of red wines has been determined.

The analysis by hyperspectral systems delivers considerable accuracy and estimation of the antioxidant activity of wines, thus proving that such technique is highly promising in determining the quantity of antioxidant activity and in the creation of classification models for sorting red wines by their antioxidant properties.